



## RAPID COVID TESTING FOR THE WORKPLACE



### ABBOT PANBIO MANUAL TESTING

### ABBOT ID NOW AUTOMATED TESTING

Antigen	Type of Test	Molecular PCR
Yes	Nasal swab administered by health professional	Yes. Swab can also be used as a throat swab.
Manual	Swab testing	Automated
15 minutes	Result Turnaround Time	10 to 12 minutes
91.1% 99.8%	Sensitivity Specificity	95% 97.9%
Between 2 and 30 degrees C	Storage	Between 15 and 30 Degrees C
No	Diagnostic machine required	Yes. Also compatible with nasal pharyngeal swab
Yes	Health Canada Emergency Use Approval	Yes
25	Tests per box	24

**Our Solutions. Your Safety.**  
**WWW.HETEK.COM**  
**1-888-432-8422**